


# iHILIC®-Fusion(P)



## Charge Modulated Amide HILIC columns

- polymer based, 5  $\mu\text{m}$
- available in PEEK
- stable at high pH
- ultra-low bleeding



| P/N          | iHILIC®-Fusion(P) HPLC Column, PEEK  |
|--------------|--|
| 130.022.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 20x2.1 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$        |
| 130.052.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 50x2.1 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$        |
| 130.102.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 100x2.1 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$       |
| 130.152.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 150x2.1 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$       |
| 130.252.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 250x2.1 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$       |
| 130.054.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 50x4.6 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$        |
| 130.104.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 100x4.6 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$       |
| 130.154.0520 | iHILIC®-Fusion(P) HILIC Column, PEEK, 150x4.6 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$       |
| 131.322.0520 | iHILIC®-Fusion(P) Guard Column, PEEK, 20x2.1 mm, 5 $\mu\text{m}$ , 200 $\text{\AA}$ , 3 pk |

**ACHROM nv : BENELUX dealer for HILICON - T + 32 9380 06 47**

**[achrom@achrom.be](mailto:achrom@achrom.be) - [www.achrom.be](http://www.achrom.be)**