

## DigiVol : Digital Programmable Syringe with a Brain



The DigiVOL is designed for precise liquid dispensing, offering exceptional reliability and accuracy in applications such as micro HPLC, solid phase extraction, filtration and membrane separation.

**Key features** of the digiVOL Programmable Digital Syringe include:

- Accurate and precise dispensing of liquids at pressure
- Programmable aspiration and dispense volumes and flow rates
- Bluetooth handle powered by an external 24-volt power supply
- Simple to operate software with Quick Run Method modes
- Rapid syringe changing with ePrep XCHANGE

**Specifications :** The digiVOL Digital syringe driver kit comes with a programmable controller, a stand, a syringe handle and a 250  $\mu$ L syringe, can be used either directly on the stand or in hand-held mode for smooth and flexible operation. The digiVOL Calibration standard kit comprises of the digiVOL digital syringe driver kit complete with the accessories that you need to make your standards quickly & easily.

[Click on DigiVol picture to read more !](#)

Description	Cat nr	Price	Now
digiVOL Digital syringe driver kit with Stand	01-08115	3500,00 €	2975,00 €
digiVOL Calibration Standard kit with Stand and 50 $\mu$ L, 250 $\mu$ L, 1.25 mL, 2.5 mL and 5 mL syringes	01-08210	4675,00 €	3975,00 €

Exclusive Promotion with Achrom valid till end of march 2024 , 21% VAT excl. and free shipping

Achrom nv – Phlipprostraat 7-9870 Belgium- [www.achrom.be](http://www.achrom.be) - [achrom@achrom.be](mailto:achrom@achrom.be) T+32 9 380 06 47