
iHILIC®-Fusion(+)



Charge Modulated Hydroxyethyl Amide HILIC Columns

- **HPLC, with enhanced retention for anions**
- **silica based, 3.5, and 5 µm**
- **available in PEEK and stainless steel (SS)**
- **high efficiency and ultra-low bleeding**

P/N	iHILIC®-Fusion(+) HPLC Column, PEEK
100.022.0310	iHILIC®-Fusion(+) HILIC Column, PEEK, 20x2.1 mm, 3.5µm, 100Å
100.052.0310	iHILIC®-Fusion(+) HILIC Column, PEEK, 50x2.1 mm, 3.5µm, 100Å
100.102.0310	iHILIC®-Fusion(+) HILIC Column, PEEK, 100x2.1 mm, 3.5µm, 100Å
100.152.0310	iHILIC®-Fusion(+) HILIC Column, PEEK, 150x2.1 mm, 3.5µm, 100Å
100.054.0310	iHILIC®-Fusion(+) HILIC Column, PEEK, 50x4.6 mm, 3.5µm, 100Å
100.104.0310	iHILIC®-Fusion(+) HILIC Column, PEEK, 100x4.6 mm, 3.5µm, 100Å
100.154.0310	iHILIC®-Fusion(+) HILIC Column, PEEK, 150x4.6 mm, 3.5µm, 100Å
100.022.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 20x2.1 mm, 5µm, 100Å
100.052.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 50x2.1 mm, 5µm, 100Å
100.102.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 100x2.1 mm, 5µm, 100Å

P/N	iHILIC®-Fusion(+) HPLC Column, PEEK
100.152.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 150x2.1 mm, 5µm, 100Å
100.252.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 250x2.1 mm, 5µm, 100Å
100.054.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 50x4.6 mm, 5µm, 100Å
100.104.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 100x4.6 mm, 5µm, 100Å
100.154.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 150x4.6 mm, 5µm, 100Å
100.254.0510	iHILIC®-Fusion(+) HILIC Column, PEEK, 250x4.6 mm, 5µm, 100Å
101.322.0510	iHILIC®-Fusion(+) Guard Column, PEEK, 20x2.1 mm, 5µm, 100Å, 3 pk

ACHROM nv : BENELUX dealer for HILICON - T + 32 9380 06 47

achrom@achrom.be - www.achrom.be